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## SALIVA (Covid19) Protein Detection Kit

EN

*For professional and in vitro diagnostic use only.*

### Intended Use

“SALIVA (Covid19) Protein Detection Kit” is a chromatographic test for the qualitative detection of SARS-CoV-2 in Saliva samples.

It provides an aid in the diagnosis of infection with SARS-CoV-2.

Positive results indicate the presence of antigen, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses.

Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results must be combined with clinical observations, patient history and epidemiological information, and confirmed with a molecular assay, if necessary for patient management.

### Principle

The kit works with the principle of antigen-antibody reaction and immunochromatography. The device includes: monoclonal antibody 1 against sars-cov-2 marker protein labeled by hypersensitive colloidal gold; monoclonal antibody 2 against sars-cov-2 marker protein is coated at T line of reaction zone; goat anti chicken antibody is coated at position of quality control line (C).

In the process of the test, when the level of coronavirus-19 in the sample reaches or exceeds the detection threshold, the antigen of coronavirus-19 in the sample binds with the monoclonal antibody 1 precoated on the gold pad. The conjugates migrate upward through capillary effect, and then bind to the coated mcab-2 at the T-line. If there is no coronavirus-19 in the sample, there will be no purple red band on the T line. No matter whether there is coronavirus-19 in the sample, purple red band will appear at the position of quality control line (C). The purplish red band of quality control line (C) can be used as the standard judge whether the sample is sufficient and whether the chromatographic process is normal or not. It can also be used as the internal control standards of reagents.

### Warnings and Precautions

- ✓ For in vitro diagnostic use only
- ✓ For healthcare professionals and professionals at point of care sites
- ✓ Do not use after the expiration date
- ✓ Do not freeze. Avoid excessive temperature and humidity in the experimental environment. The reaction temperature should be 15°C-30°C and the humidity should be below 70%
- ✓ The package bag contains desiccant and it should not be taking orally
- ✓ Please read all the information in this leaflet before performing the test
- ✓ The test cassette should remain in the sealed pouch until use
- ✓ All specimens should be considered potentially hazardous and handled in the same manner as an

infectious agent

- ✓ Do not eat, drink, smoke or chew gum etc. within 30 minutes before giving a saliva sample
- ✓ Wash with water if stabilizing liquid contracts with eyes or skin
- ✓ Tightly cap the tube after collection, ensure on leakage
- ✓ The used test cassette should be disposed of as hospital waste material

## Composition

Each kit contains a test cassette and a bag of desiccant, a set of saliva collector including a saliva funnel and a collection tube containing 1ml diluent and a dropper.

The test cassette consists of gold label pad, sample pad, nitrocellulose membrane, absorbent paper, PVC board and plastic card.

## Kit components


- ✓ 1 Disposable test card
- ✓ 1 saliva funnel
- ✓ 1 collection tube containing 1 ml diluent
- ✓ 1 Dropper

## Storage and Stability

- ✓ Store as packaged in the sealed pouch at the temperature (2-30°C). The kit is stable within the expiration date printed on the labeling
- ✓ Once open the pouch, the test should be used within 30 minutes after unsealing. Prolonged exposure to hot and humid environment will cause product deterioration
- ✓ The LOT and the expiration date were printed on the labeling

## Test Procedure

Allow the test device and specimens to equilibrate to temperature (15-30°C) prior to testing.

1. Don't eat, drink, brush your teeth or use mouth wash, etc. before 30 minute using the test.
2. Open the aluminum foil bag of the test cassette, take out the test cassette and mark the inspected person or sample number on the cassette. Use within 30 minutes, especially at room temperature above 30 °C or high humidity, as soon as possible.
3. Put the kit on a clear platform, open the cover of the diluents tube, screw on the saliva funnel, make a <Kuuua> sound in the throat to clear Saliva to the deep throat, and collect the deep throat saliva to 2ml

(there is already a 1ml solution in the diluents tube); screw off the saliva funnel, cover it, turn upside down and mix well for 6-10 times, then screw off the cover, let stand for 10 seconds, suck the liquid with a dropper, drop 3-4 drops (no bubble and with 3 seconds interval between each drop) into the sample hole, and start to count for 10-15 minutes.

4. Wait for the purple red stripe to appear. The test results should be read within 10-15 minutes, after more than 15 minutes, the reading results are invalid.

## Interpretation of results:

**Positive(+):** Purple-red bands should be appear in C and T areas

**Negative(-):** Purple-red band should be appear only in C area

**Invalid:** No matter whether there is purple red band at T line position, if there is no purple red band at quality control line (C), it indicates that the operation procedure is incorrect or the test kit has been invalid.

## Sample Transport and Storage

Samples should be used as soon as possible after collection and should not be stored for a long time at room temperature. If the sample can not be detected in time, the sample can be stored for 48 hours at 2 °C – 8 °C. For long term storage, sample should be frozen at – 20 °C, avoiding repeated freezing and thawing.

## Limitations

- ✓ “SALIVA (Covid19) Protein Detection Kit” is a qualitative and preliminary screening methods, which can not determine the content of coronavirus-19 in saliva. This kit provides only one preliminary analysis result, and a second analysis method must be used to determine the result.
- ✓ The possible causes of negative results of this product include:
  - There is no virus in the sample, or the virus content in the sample is very low, which is lower than the critical detection concentration of reagent, and the low concentration sample cannot be detected.
  - The virus may have been inactivated in the sample. Although the nucleic acid fragment may still exist, the virus antigen may have been destroyed and inactivated.
  - Incorrect operation or other factors that may affect the detection, such as improper transportation and storage, leading to reagent failure.
- ✓ When the epidemic degree of the disease decrease, the positive predictive value decreases, so the interpretation of positive results of low-risk population should be cautious.

## Performance characteristics

### Cross Reactivity (analytical specificity)

“SALIVA (Covid19) Protein Detection Kit” has no cross reactivity with endemic human coronavirus HKU1, NL63, 229E, adenovirus, influenza A, influenza B, measles virus, rotavirus, norovirus, mumps virus, human cytomegalovirus, varicella zoster virus, respiratory syncytial virus, Chlamydia pneumonia, mycoplasma pneumonia.

## Interference

“SALIVA (Covid19) Protein Detection Kit” results do not be interfered with the substance at the following concentration: bilirubin concentration  $\leq 250\mu\text{mol/L}$ , triglycerides concentration  $\leq 15\text{mmol/L}$ , hemoglobin concentration  $\leq 9\text{g/L}$ .

“SALIVA (Covid19) Protein Detection Kit” results do not be influenced by the following substance:  $\alpha$ -interferon, zanamivir, ribavirin, oseltamivir, paramivir, lopinavir, ritonavir, abidol, levofloxacin, azithromycin, ceftriaxone, meropenem, tobramycin, histamine, hydrochloride, phenylephrine, oxymetazoline, sodium chloride (containing Preservatives), beclomethasone, dexamethasone, flunisolide, triamcinolone, budesonide, mometasone and fluticasone.

## Clinical Performance

Method	RT-PCR		Total	
	Results	Positive		Negative
SALIVA (Covid19) Protein Detection Kit	Positive	51	0	51
	Negative	1	50	51
Total		52	50	102

- ✓ Sensitivity : 98.07%
- ✓ Specificity : 100.00%
- ✓ Accuracy : 99.01%

## Symbols used

	Do not reuse		For in vitro diagnostic use only
	Store between 2-30°C		Consult instructions for use
	Caution		Lot number
	Use by		Contains sufficient for <n> tests
	Keep away from sunlight		Keep dry
	Manufacturer		Do not use if package is damaged

## SALIVA (Covid19) Protein Detection Kit

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