



GRAPADORA DE PIEL DESECHABLE



DESCRIPCIÓN DEL PRODUCTO

REF: KHPF-35A

- Las grapas son compatibles con el tejido de la piel, **sin reacción de rechazo**.
- Soportan grapas de diferentes dimensiones, se aplican para grapar diferentes heridas de la piel.
- Mucho más rápido que la sutura a mano, por lo que alivia la carga de trabajo de la operación, acorta el tiempo de la misma y disminuye la aparición de complicaciones.
- Evitan dañar el tejido de la piel, lo que beneficia la curación de las heridas. Grapa de manera suave y precisa, la cicatriz de curación es pequeña y limpia, cuando se quitan las grapas, no hay sensación de dolor. La **cabeza angulada** de la grapadora proporciona una **visión clara para asegurar un grapado** preciso.

ESPECIFICACIONES DEL PRODUCTO

- Altura de la grapadora cerrada: 7.0(A) x 4.5(A)mm
- Material de la grapadora: Acero inoxidable médico

Unit: mm

Model	Spec.	Type	Staples								Number of Staples	Optional (Skin Staple Remover)
			Diameter	Tolerance	Dimensions after forming (L1 x H1)	Tolerance	Dimensions before forming (L x h)					
							L	Tolerance	H	Tolerance		
KHPF	10, 15, 25, 35, 45	W	0.6	±0.1	7.3 x 4.0	±0.6	13.1	±0.2	3.5	±0.1	10, 15, 25, 35, 45	KHPF-CA
	A	0.6	7.0 x 4.5		13.1		3.5		KHPF-CB			
	B	0.6	5.5 x 3.5		11.3		2.7					