

Disposable Monopolar Instrument



- Convenient, single use instrument, eliminating processing and sterilization costs.
- Extended insulation to prevent stray currents.
- Ergonomic handle for surgeon comfort and control.
- 360 degree finger wheel rotation for simple, precise jaw alignment.



Product Code	Product Discription	Size
MMS205330	Disposable Monopolar Curved Scissors 17mm Blade	5mm*33cm
MGP305330	Disposable Monopolar Grasper Forceps (Ratchet)	5mm*33cm
MCH405330	Disposable Monopolar Clinch Forceps (Ratchet)	5mm*33cm
MBK505330	Disposable Monopolar Babcock Forceps (Ratchet)	5mm*33cm
MUG605330	Disposable Monopolar Uterus Grasping Forceps (Ratchet)	5mm*33cm
MDS705330	Disposable Monopolar Dissector Forceps (NO Ratchet)	5mm*33cm
MDM805330	Disposable Monopolar Mixer Forceps (NO Ratchet)	5mm*33cm
MAI905330	Disposable Monopolar Atraumatic Intestinal Forceps (Ratchet)	5mm*33cm

SOLD 10 UNITS PER BOX